



**Shape Memory  
Aligner®**

# **World's 1<sup>st</sup> Direct Printed Aligner Shape Memory Aligner®**

World's first patented technology Shape Memory



[www.graphy3d.com](http://www.graphy3d.com)  
[www.itgraphy.com](http://www.itgraphy.com)


**Graphy**

3D Print the World with Graphy's Solutions

Eng

# About Graphy Inc.

Korea's leading manufacturer of 3D printer materials



## Graphy Inc.

Is a company that develops and manufactures new materials (photocurable resin) that specializes in 3D printing-based technology, a major element of the 4th industrial revolution, and has optimized knowledge of new material technology and medical convergence technology. With the successful development of Tera-Harz, we are growing into a global company through partnerships with the world's leading companies in the 3D printer material business field.

Graphy is not going to stop at manufacturing and developing materials, but aims to apply 3D printer technology to real life by supplying various 3D printer application solutions and medical solutions to the market.

## Conference & Exhibition



# Shape Memory Aligner® Specialists Around the World

What world-class scholars and speakers are saying about Shape Memory Aligner®



**Prof. Ravindra Nanda Professor Emeritus**

Department of Orthodontics University of Connecticut  
Adjunct Professor The Forsyth Institute  
Editor-in-Chief Progress in Orthodontics

"What made the big difference in orthodontics especially with the wires and bracket was when we introduce shape memory wires, wires which change their characteristic with the heat as well as when you're going to put up in the mouth then force is going to remain pretty stable over a long period of time. And this is what now Graphy aligner material is going to do."



**Prof. Ki Beom Kim, DDS, MSD, PhD**

Dr. Lysle Johnston Endowed Chair in Orthodontics  
Professor and Program Director  
Department of Orthodontics  
Center for Advanced Dental Education  
Saint Louis University

"I think the Graphy material is an innovative material that can fundamentally overcome the disadvantages and limitations of existing thermoplastic materials, as well as the simplicity of device manufacturing."



**Dr. Werner Schupp**

Editor in Chief "Journal of Aligner Orthodontics" Quintessence Publishing  
Editor for "Manuelle Medizin und Kieferorthopädie" (Springer Verlag)  
Advisory Board "Kieferorthopädie" Quintessenz Verlag

"Advantages of 3D aligner printing: Tera Harz Clear is more flexible and has a larger elastic range. Perfect fit to tooth, almost no gap"



**Dr. Kenji Ojima**

D.D.S M.D.Sc  
Invisalign faculty  
Adjunct Professor University of Torino  
President of Japan Academy of Aligner Orthodontics

"This material is a very new future. There's a better way for new generation aligners on orthodontics."



**Dr. Björn Ludwig**

Assistant Professor at the University of Homburg/Saar,  
Department of Orthodontics  
Editor in chief of the Quintessenz publication "Kieferorthopädie"  
Co-editor of the Journal of clinical Orthodontics

"I had to print out the aligners the same day, and it worked, because I can print it out in my office. It gives me control, it gives me the power."



**Dr. Yong-Min Jo**

Since 2020 Shareholder and Chief Business  
Development & Innovation Office at Scheu Group:  
Scheu Dental, CA Digital, Smile Dental Companies

"The innovative advantage of the new aligner material is the continuous power transmission thanks to its memory effect. Due to this effect, the aligner permanently returns to its original shape, which has a positive effect on the power transmission and rigidity of the aligner."



**Dr. Nearchos C. Panayi**

Scientific coll. European University Cyprus  
Visiting Research Scientist, Clinic of Orthodontics, University of Zurich,  
Author "Design It Yourself Orthodontics" book, Quintessence  
Associate Editor "Progress in Orthodontics"

"Preliminary results show excellent mechanical properties. It is clear to see that aligner printing is here to stay. The river cannot reverse. Nevertheless, in order to have a consistence and excellent result, Graphy has created a protocol which needs to be followed beginning from the appropriate aligner design, printing, and most important, UV curing."



**Dr. Simon Graf**

Smile AG  
Inventor and pioneer in CAD/CAM  
procedures for metal and acrylic 3D printed  
orthodontic appliances

"The two ways of aligners!  
Now direct printed instead of vacuumforming."



**Dr. Akim Benattia**

Dr. Akim Benattia has served as a former assistant at universities in France and currently holds positions as the General Secretary of the French Aligner Society, as well as memberships in The European Aligner Society, the Japan Aligner Academy, and the Deutsche Gesellschaft für Aligner Orthodontie. Since 2015, he has been an international speaker, presenting in countries such as France, Tunisia, Italy, Spain, Morocco, Taiwan, Ukraine, Greece, Dubai, Latvia, and Austria.

Advancements in Graphy materials significantly enhance the potential of Clear Aligner Therapy thanks to its Shape Memory properties



**Dr. Choon Gwack**

PNU School of Dentistry DDS PhD  
Vice President of the KAO Busan Gyeongnam-Ulsan Division  
Vice President of the KSDO

Orthodontic practice director of BarunE Orthodontic Clinic  
"Made from shape memory polymers, Shape Memory Aligner®s are moderately elastic at intraoral temperatures and soften at higher temperatures. They have the advantage of being manufactured without blockouts for easy insertion, so you can achieve good, predictable alignment with minimal attachments."



**Dr. Young Jin Jeon**

Graduate of PNU School of Dentistry, DDS PhD  
Invisalign Speaker for Asia-Pacific (APAC)  
Assistant Professor of Orthodontics, Case Western Reserve University  
Orthodontic Practice Director, Erumi Orthodontic Clinic

Just as Shape Memory Wire enhanced bracket orthodontics, Shape Memory Direct Aligner will revolutionize clear aligner treatment.

# What is Shape Memory Aligner®? World's 1<sup>st</sup> 3D Printed Shape Memory Aligner®

## Shape Memory Aligner® is

more than just a clear aligner.

It's not any typical clear aligner, and it's not just a 3D print.


It is the world's first patented product to feature Shape Memory.

A state-of-the-art aligner recognized by orthodontic masters around the world.


The Shape Memory Aligner® is a new paradigm in orthodontics, far beyond brackets and other traditional methods. Shape Memory Aligner® is an aligner that uses digital information from a patient's mouth to print directly to a 3D printer, bypassing the step of creating a tooth model.

Shape Memory Aligner® is a medical device that reduces production time, is suitable for biomechanical movements such as adjusting the thickness of each tooth, which is not possible with conventional thermoforming manufacturing methods, and is harmless to the human body.

## Shape Memory Aligner® mechanical property



Shape memory effect and shape recovery ratio over time of TC-85  
04




Reference: Lee, S. Y., Kim, H., Kim, H. J., Chung, C. J., Choi, Y. J., Kim, S. J., & Cha, J. Y. (2022). Thermo-mechanical properties of 3D printed photocurable shape memory resin for clear aligners. Scientific Reports, 12(1), 1-10.

## Shape Memory Aligner®


## Thermoformed Aligner

## Adjustable strain and resilience




※ Based on self-test of Graphy

## Shape driven biomechanics




However, with the Shape Memory Aligner®, the released forces keep coming back over and over again by the effect of Shape Memory.


Reference: Koenig, N., Kim, KB. D.D.S., M.S.D., Ph.D. (2022). Comparison of dimensional accuracy between direct-printed and thermoformed aligners. The Korean Journal of Orthodontics, Article In Press. Published online



※ Based on self-test of Graphy




This is a displacement driven system and not a force driven system, but Shape driven biomechanics is ideal.



The deformation of the arch form of the Shape Memory Aligner® does not cause deformation of the entire tooth span.

In addition to the risk of fracture when 3D printing aligners with traditional materials, aligners without shape memory limit the orthodontic ability to return to its original shape when it is deformed



When the center red circle loses strength for the limit of the material's stress, the orthodontic force weakens.

※ Based on self-test of Graphy

## Shape Memory Aligner® patents and medical device safety certifications

### Direct printing allows for strategic thickness control of orthodontic treatment

The thickness of the Shape Memory Aligner® is adjustable part by part, allowing for the fabrication of customized treatment plans for each patient. On the other hand, conventional clear aligners use thermoformed sheet paper, which is easily deformed. The compression thermoforming makes it difficult to produce aligners of the intended thickness, and it is impossible to adjust the thickness of each tooth to provide orthodontic force to the tooth, resulting in cumulative errors due to multiple steps in the process.

Shape Memory Aligner®

Easy brushing

Healthy teeth and gums, deodorizing bad breath


Maximize the effectiveness of orthodontic treatment

Less stress on the jaw joint

## Minimized inconvenience when putting on and taking off, hygienic management through heat disinfection

Minimizing block out for intraoral wearing and removal improves orthodontic power and makes patient more comfortable. The Shape Memory Aligner® can be heat disinfected by placing it in hot water, so it's always in pristine condition. It's also made of body-friendly materials.

## Light & Continuous force



Aligners made of conventional thermoplastic sheets generate an initial heavy force that strains the teeth, causing patient pain, and then rapidly drops below the appropriate force, resulting in an inability to maintain the tooth movement force and achieve the intended orthodontic result.

Due to its innovative technology and material properties, the Shape Memory Aligner® allows for light & continuous force to achieve the intended alignment results. Shape Memory Aligner® induces light and sustained force. Light & Continuous Force induces smooth, continuous tooth movement while activating cells.

Reference : Contemporary Orthodontics. William R. Proffit et al. 6 edition p.263-266

Application No.	Patent No.	Title of Patent
10-2018-0099277	Patent 10-2067533	RADIATION CURABLE POLYMER COMPOSITION FOR 3D PRINTER
10-2018-0123132	Patent 10-2111789	PHOTOCURABLE COMPOSITION FOR 3D PRINTER FOR PRODUCING TRANSPARENT ORTHODONTIC DEVICE
10-2020-0033430	Patent 10-2281789	Patient-specific orthodontic device using 3D printing and orthodontic method using the same
10-2020-0172326	Patent 10-2244439	Transparent aligner with increased correction power
10-2021-0052197	Patent 10-2273407	Post-curing method of 3D printer output and transparent orthodontic device manufactured by the method

World's first "FDA 510k Class II", "CE Class II", "MFDS Class 2", "PMDA Class II" Cleared



## Shape Memory Aligner® Related articles




Shape Memory Aligner® about 50 relevant SCI-quality international journal articles  
Collaborations with overseas companies and research institutes and schools are expanding and are now academically recognized.

# Orthodontics platform

## DAD In-house Producuion

- Graphy's 'Direct Aligner Designer (DAD)' is orthodontic simulation software that allows treatment plan establishment, setup, and design work to be performed in one software.
- The automatic separation and alignment function can provide orthodontic simulation results for patients with malocclusion within 5 minutes.



## GraphyTeams Interactive platform for Shape Memory Aligner

- Graphy provides a custom production service where users do not produce shape memory aligners themselves, but instead request production and have it delivered.
- 'GraphyTeams' allows you to intuitively check reception details, progress status for each patient, delivery status, etc. through a cloud-type database.



www.graphyteams.com



## THC 3rd Generation Nitrogen Curing Machine


# Tera Harz Cure

World's first UV curing device to achieve 100% polymerization conversion of materials


- 280,000Mj/cm2,1,000MW/cm2(5 minute-curing based), the best material intensity and shade with high-intensity UV energy
- Ability to store 5 curing conditions per material
- The optimal light level for each printed materials(Level 1~5)
- 360° UV irradiation and optimal LED arrangement for best light uniformity
- Safety cooling system based on LED temperature
- Maximum curing size(180mm,360 Turntable)



[ Nitrogen Generation ]



[Minimize Oxygen in the Chamber (artificial image)]



### Benefits of Nitrogen Curing

- Able to obtain the best mechanical properties and optimal shading
- Minimize oxygen to improve surface quality for C&B, denture fabrications, etc.
- Perfect final metal crown alignment by improving casting surface quality
- Nitrogen curing is mandatory for dental and medical field products (for surface tacky and residual polymer removal)
- Improved economic efficiency and convenience compared to existing nitrogen curing machines connected to a nitrogen tank.

Properties	Description
Display	7.9" TFT Touch LCD
LED Wavelength	405nm
LED Power	200W
UV Energy Density	280,000 mJ/cm²
Irradiance of UV (5 minute-curing based)	1,000 mW/cm²
Curing Chamber	Ø180 × 650mm
Dimension (weight)	275 × 310 × 310mm (8.5kg)

# Tera Harz Smart Robot

The most advanced integrated device for  
Shape Memory Aligner manufacturing.



## Advantages of Tera Harz Smart Robot

This system is more than just a system.  
The excellence of the world's first Shape Memory Aligner is more elaborated, possible by introducing  
a state-of-the-art robot-system for greater excellence and convenience.  
Shape Memory Aligner is made more efficient with the solution of Graphy's new Tera Harz Smart Robot.

- 01

**Competitive excellency in technology**  
Deliver consistent quality of product, 24hr operation time and optimized automated reprogramming.
- 02

**Improvement outcomes**  
Improve the quality of the outcome through sophisticated operation.
- 03

**Economic benefits**  
Reduce consumables, accelerate the treatment and enhance productivity.
- 04

**Time efficiency**  
Same day and convenient treatment.

Product Name	Tera Harz Smart Robot (THSR)
Model Name	R2K2U
Weight(Kg)	20
Application	In-house Aligner Manufacturing (lab/clinic)
Voltage	220V(110V Compatible Module Provided)
Usage	Production Line
After Warranty Service	Video Technical Support
Type	6-axis Vertical Multi-joint
Machine type	Aligner Manufacturing Robot Arm 6 Axis
Aligner Material	Tera Harz Clear
Warranty	1 year

# World's 1<sup>st</sup> Shape Memory Aligner®









## Shape Memory Aligner® Package

Treatment Option	SMA Lite II	SMA Standard	SMA Advanced
Number of Steps	Up to 16	Up to 30	More than 30
Treatment Period	4~6 Months	7~10 Months	3 Years
Warranty Period	12 Months	2 Years	3 Years
Additional Extension Aligner	N/A	Max. 15 Steps	No limit within warranty period
Scope of Treatment	Normal Orthodontic Treatment	Normal Orthodontic Treatment	Extraction Case, Complex treatment
Treatment Example	<div>■ <b>Anterior</b></div> <div>A. Crowding: Simple ~ Moderate</div> <div>B. Spacing case</div> <div>C. In/Extrusion: under 1-2mm</div> <div>D. Rotation (under 15 °)</div> <div>E. Relapse case</div> <div>■ <b>Posterior</b></div> <div>A. In/Extrusion : less than 1mm</div> <div>B. Distalization</div>	<div>■ <b>Anterior</b></div> <div>A. Severe Crowding</div> <div>B. 4 anteriors Rotation: under 30°</div> <div>C. In/Extrusion : over 2mm (elastic)</div> <div>■ <b>Posterior</b></div> <div>A. Canine – Posterios rotation : under 15°</div> <div>B. In/Extrusion in 2mm (Attachment/Elastic)</div> <div>C. Mesio/Distalization</div>	<div>■ <b>Anterior</b></div> <div>A. Extraction case - anteriors/ premolar</div> <div>B. Easing skeletal malocclusion</div> <div>C. Rotation : over 30°</div> <div>D. Posterior expansion</div> <div>E. In/Extrusion : over 2mm(Att/ Elastic)</div>

## Shape Memory Aligner® Case







Initial		Final
	<b>Lite</b>  Age: 33Y Gender: Female Tx Period: 10 months No. of steps: 37	
		
		

Courtesy of **Graphy**







Initial		Final
	<b>Lite</b>  Age: 16Y Gender: Male Tx Period: 6 months No. of steps: 18	
		
		

Courtesy of **Graphy**

## Shape Memory Aligner® Case

Initial		Final
	<b>Standard</b>  Age: 55Y Gender: Male Tx Period: 7 months No. of steps : 27	
		
		

Courtesy of **Dr. C Gwack**

Initial		Progress
	<b>Advanced</b>  Age: 25Y Gender: Male Tx Period: 1 year No. of steps: 49  Extraction Surgery-first Approach	
		
		

Courtesy of **Dr. ES Kim**

## Shape Memory Aligner® Case | Advanced

### Moderate to severe crowding case

Initial



Final




Achieving Excellent Outcomes in Narrow Dental Arch Cases  
Without Extractions

Courtesy of Dr. Kenji Ojima, Japan

### Asymmetric premolar extraction

Initial



Final




Outstanding Orthodontic  
Results in Asymmetry  
Case by  
Asymmetric Extraction

Courtesy of Dr. Kenji Ojima, Japan


## Shape Memory Aligner® Case | Advanced

### MSE [Maxillary skeletal expander] combined case


Initial




Progress




Final



Courtesy of Dr. Kenji Ojima, Japan



Advances in orthodontic clear aligner materials.  
Bioactive Materials, 22, 384-403.



Courtesy of **Dr. Simon Graf, Swiss**



### In-House 3D-Printed Shape Memory Aligners for Retreatment after Fixed Retainer Failure

MARTIN SAUNDERS, *Spence, CO, PhD*  
 JENNIFER J. DAVIS, *Spence, PhD*  
 SARAH K. MARRAS, *Spence, MS, PhD*  
 Tel.: 719/334-7373; Fax: 719/334-7370

**R**ecent advances in materials science and computer-aided design and manufacturing (CAD/CAM) technology in combination with the growing number of adult patients seeking inconspicuous orthodontic treatment,<sup>1</sup> have led to the widespread adoption of clear aligners as a comprehensive alternative to fixed appliances. Increasingly, the digital design and fabrication of clear aligners are occurring "in-house"—within a private practice.<sup>2</sup> The introduction of the first resin specifically intended for clear aligners, one also known to be 3D-printed with greater dimensional precision than thermofused aligners,<sup>3</sup> while incorporating features such as planned variations in thickness across the aligner body and thermomechanical shock memory,<sup>4,5,6</sup>



© 2014 John Wiley & Sons, Ltd. *J. Clin. Neurophysiol.* 35, 1–11 (2014)



## Shape memory aligners: A new dimension in aligner orthodontics

**NOTES:** Drawing and structure, and engineering design, copyright © 2000 by Autodesk, Inc. All rights reserved. Autodesk reserves the right to alter product offerings and specifications at any time without notice, and is not responsible for typographical or graphical errors that appear in this document.

[illegible][illegible]

**Abstract**

combination of 2D lingual appliance and in-house aligner: An anterior cross-bite and TMD case report after 1-year follow-up



**T. G. Campbell, author of  
"Fingerprinting  
the Brain: The Neuroscience  
and Rehabilitation of Speech  
and Language"**

*Dr. Campbell is a professor of  
Psychiatry and Behavioral  
Sciences at the University of  
California, San Francisco.*

**Journal of Speech, Language  
and Hearing Disorders**  
Volume 18 Number 1  
Spring 1983

**0363-2764/83/0018-0000\$01.00/0**

**INTRODUCTION**

Intercine ester lase is defined as a multicatalytic, sodium-free, the largest performing of the microbial ester lase. It is substantially to the microbial ester lase [1]. Several cases of intercine lase having a high activity of a weak acid and also non-mutated lase [2]. Previous researches have shown that intercine lase is a weak acid and also non-mutated lase [2]. Previous researches have shown that intercine lase is a weak acid and also non-mutated lase [2]. Previous researches have shown that intercine lase is a weak acid and also non-mutated lase [2].



Fig. 1 32-year-old female patient with severe protrusion of lower left central incisor, lateral incisor, and canine due to failure of fixed lingual retainers before retreatment.



Fig. 4 After 24 weeks, including nine steps of treatment and three steps of refinement.




Figs 18a-l Comparison (a-c) before and (d-f) after orthodontic treatment with in-office direct Graphy aligners, (g to i) After restoration with composite on the maxillary and mandibular anterior teeth and 1 week after gingivectomy (Dr Wolfgang Boisseree, Cologne). (j to l) Smile before and after interdisciplinary treatment.



**Figure 2: Pre-treatment X-rays.**

Figure 9: Post-treatment X-rays



[itgraphy.com](http://itgraphy.com)



[graphy3d.com](http://graphy3d.com)

Graphy Inc. 6th Fl. AceGasamFORHU, 225, Gasan digital 1-ro, Geumcheon-gu, Seoul, 08501 South Korea

**TEL** +82-2-864-3056 | **FAX** +82-2-864-3057 | **E-mail** [info@itgraphy.com](mailto:info@itgraphy.com) | **Web** [itgraphy.com](http://itgraphy.com), [graphy3d.com](http://graphy3d.com)

Revised Version Published September 2024